

Spectrophotometry

TASK: Determine by spectrophotometer "JanyWay" the concentration of riboflavin vitamin B₂ and appreciate sensitivity calculating molar absorbtion coefficient:

Nr.1, Nr.2, Nr.3, Nr.4, Nr.5, Nr.6., Nr.7, Nr.8, Nr.9, Nr.10, Nr.11,Nr.12., Nr.13, Nr.14.

Calibration graph preparation history.

1. Standard solution of riboflavin $C_{\text{riboflavin}}=0.04 \text{ mg/mL}$ have been prepared in 8 test tubes 10 mL standard solution with distilled water.

2. Calculated concentrations in 8 tubes are fix results in table: $C_{B_2} = \frac{C_{\text{Riboflavin}} \cdot V_{B_2}}{10\text{mL}}$.

3. Choose on spectrophotometer "JanyWay" wave length $\lambda = 445 \text{ nm}$ and measure absorption $A_x = \log(I_0/I)$ of chosen **samples** Nr. __. Use calibration graph from table and read unknown

concentration of the **sample** C_x . Table of results

<i>Solution Nr.</i>	<i>blank</i>	<i>1</i>	<i>2</i>	<i>3</i>	<i>4</i>	<i>5</i>	<i>6</i>	<i>7</i>	<i>8</i>
$V_{B_2} \text{ (mL)}$ $C_{B_2}=0.04$	0.00	0.50	1.00	1.50	2.00	2.50	3.00	3.50	4.00
mL $V_{\text{ūdens}}$	10.00	9.50	9.00	8.50	8.00	7.50	7.00	6.50	6.00
Absorption $A = \log(I_0/I)$	0.000	0.063	0.121	0.184	0.247	0.302	0.368	0.426	0.482
$C_{B_2}, \text{mg/mL}$	0,000	0,002	0,004	0,006	0,008	0,010	0,012	0,014	0,016

Calculate the molarity of **B₂** solution $C_M = C_x / M_{B_2}$ and molar absorption factor $a \text{ M}^{-1}\text{cm}^{-1}$, if glass cell thickness size is measured. Use the ruler and check the glass cell thickness! $l = 1 \dots \dots \text{ cm}$ or no?

$$C_M = C_x / M_{B_2} = \dots \dots \dots \text{g/L} / \dots \dots \dots \text{g/mol} = \dots \dots \dots \text{ M}^{-1}$$

$$a_{445} = A_x / C_M \cdot l = \frac{A_x}{C_M \cdot l}, \text{ M}^{-1}\text{cm}^{-1} = \dots \dots \dots = \dots \dots \dots \text{ M}^{-1}\text{cm}^{-1}$$

Concentration of **sample** : $C_x = \dots \dots \dots \text{ mg/mL} = \text{g/L}$

Graph for calibration line $A = a \cdot C \cdot l$ drawn through zero **0** value . $A = C = 0$

On spectrophotometer "JanyWay" measure the **sample** for analyze absorption A_x and calculate its concentration C_x using the graph (**mg/mL**).


Give the calculated results for riboflavin content **mg%**, which shows

$$\text{mg}\% = \text{mg} / 100\text{mL} \cdot 100\% = \dots \dots \dots \text{mg}\%$$

The molar mass of oxidized form riboflavin (B₂ vitamin) is $M_{B_2} = 454.35 \text{ g/mol}$.

$$\uparrow A = a \cdot C \cdot \ell = \log(I_0/I)$$

Place for reading data from calibration graph A - C_{B_2}



Use calibration graph straight line through zero point for each sample concentration in units of mass concentration
 $\gamma = C_{\text{riboflavins}}, \text{ mg/mL=g/L!}$

1) Determine measured samples mass concentration, using calibration straight line!

3. table. Measured sample solutions concentrations

Sample Nr.	1	2	3	4	5
Absorption A_x	0,147	0,452	0,086	0,228	0,280
Concentration γ , mg/mL	0,0046	0,014	0,00275	0,0075	0,00925

Calibration graph straight line $A = a \cdot C \cdot \ell$ mathematically through zero point, as both are $A=C=0$.

On "JanyWay" spectrophotometer measure **sample** absorption A_x and use the graph for its concentration C_x (mg/mL) evaluation.

Converting riboflavin content to milligram percents mg%, what shows amount on 100 mL samples
 $\text{mg}\% = \text{mg}/100\text{mL} \cdot 100\% = \dots\dots\dots \text{mg}\%$

Mol mass oxidized riboflavin B₂ vitamin form is $M_{B_2} = 454.35$ g/mol.

Calculate molarity of **B₂** vitamin solution $C_M = C_x / M_{B2}$ and molar absorption coefficient $a \text{ M}^{-1}\text{cm}^{-1}$

¹, if cuvette thickness measured. With ruler check the thickness! $l = 1 \text{ cm} !?$

$$C_M = C_x / M_{B2} = 0.017 \text{ g/L} / 454.35 = 3.74 \cdot 10^{-5} \text{ (mol/L) M}$$


$$a_{445} = A_x / C_M \cdot l = \frac{A_x}{C_M \cdot l}, \text{ M}^{-1}\text{cm}^{-1} = 0.53 / 3.74 \cdot 10^{-5} = 14171 \text{ M}^{-1}\text{cm}^{-1}$$

Absorption A_x give concentration: $C_M = A_x / a_{445} / l = 0.53 / 14171 / 1 = 3.74 \cdot 10^{-5} \text{ (mol/L) M}$

Figure. Riboflavin **B₂** vitamin oxidized and reduced forms specters
Molar absorption coefficients $a_{360} = 10500 \text{ M}^{-1}\text{cm}^{-1}$ at $\lambda = 360 \text{ nm}$ and
 $a_{445} = 15499 \text{ M}^{-1}\text{cm}^{-1}$ at $\lambda = 445 \text{ nm}$.


Proteins sometimes decrease $pK_a = 10.3$ **N(3)-H** proton
protolysis constant of free flavin and decrease molar
absorption coefficient $a_{445} = 9200 \text{ M}^{-1}\text{cm}^{-1}$. For example in
glycolate oxidase (4).

Proteins cause bound flavin protolytic deprotonation
 $\text{N(3)-H} \Rightarrow \text{N(3)}^- + \text{H}^+$. Latin *flavus* – yellow.



FMN Flavin mono nucleotides riboflavin-5'-phosphate **B₂** vitamin

<http://aris.gusc.lv/FlavinMonoNucleoB2vitamPO4.Tgf>



FMN + 2 H oxidized form **Ox** + **2 H** = Reduced form **RedH₂ FMNH₂**

Oxidized form bind two hydrogen atoms +**2 H** turns to reduced form **RedH₂**,

Both and **RedH₂** forms are water soluble two hydrogen atoms **H** and **H** carrier
which transfer atoms together with two electrons as sum $2H^+ + 2e^-$.

Biochemical **oxidation – reduction** reactions progressing with water soluble
two electron equivalents transfer.

Biochemistry call vitamins about enzyme cofactors.

Small molecules for transfer call about cofactors. Vitamin **B₂** is cofactor for oxidoreductases.

Cofactors - vitamins will study in course of Biochemistry.

Conclusions

1. Riboflavin B₂ vitamin is 2 H carrier in ENZYME class E....
- 1b. Reduction is: hydrogen addition..... 1c. Oxidation is hydrogen removing.....
2. Light absorption for riboflavin are at wavelengthnm and.....nm
3. Ratio falling light intensity **I₀** over throughout going light **I** as logarithm
absorption measure is **A_x**=.....
4. Light absorption calculates as logarithmic expression **A_x=log(I₀/I)**=
5. Riboflavin concentration in the sample **Nr**..... exhibits absorption **A_x**=.....
6. Beer-Bugeers-Lambert's Law **A_x = a₄₄₅•C• ℓ=log(I₀/I)** according light absorption **A_x** is
proportional to vitamin B₂ concentration **C_x**=.....g/L
Riboflavin Vitamin B₂ molar concentration is **C_M**=.....mol/L.
7. Calculated molar absorption coefficient at λ=445 nm **a₄₄₅**=14171

Reaction velocity proportional on compounds concentration

Task: Observe the influence of $\text{H}_2\text{S}_2\text{O}_3$ sodium thio-sulfuric acid concentration $\text{CH}_2\text{S}_2\text{O}_3$ to rate of



$$\vec{v} = \frac{\Delta[\text{S}]}{\Delta t} \quad , \quad (2)$$

$$\vec{v} = k \cdot \text{CH}_2\text{S}_2\text{O}_3 \quad (3)$$

$\text{H}_2\text{S}_2\text{O}_3$ arises due to reaction: $\text{Na}_2\text{SN}_2\text{O}_3 + \text{H}_2\text{SO}_4 \Rightarrow \text{H}_2\text{S}_2\text{O}_3 + \text{Na}_2\text{SO}_4$, (4)

where $\text{H}_2\text{S}_2\text{O}_3$ is weak acid electrolyte $\alpha \rightarrow 0$ that decays according reaction (1).

Description of work: Fill in three identical test tubes (§ 1, § 2, § 3, § 4, § 5) 0.5 N $\text{Na}_2\text{S}_2\text{O}_3$ sodium thio-sulfate solutions plus distilled water (pure H_2O) with in different concentrations of $\text{Na}_2\text{S}_2\text{O}_3$;

2.) In other five test tubes 2 mL 1.5 N H_2SO_4 sulfuric acid.

In § 1 test tube mix $\text{Na}_2\text{S}_2\text{O}_3$ with 2 mL 1.5 N H_2SO_4 and switch on timer!

I ----- II ----- Table of DATA $\text{Na}_2\text{S}_2\text{O}_3 + \text{H}_2\text{O} + \text{H}_2\text{SO}_4$

Salt Nr.	$\text{Na}_2\text{S}_2\text{O}_3$ 0,5 mol/L volume, mL	H_2O mL	Create solution $\text{Na}_2\text{S}_2\text{O}_3$ concentration, mol/L	Time t, seconds	velocity, $\vec{v} = 1/t, \text{ s}^{-1}$	Time for concentration, $t_C = \tau^{1/2} \ln(C^0/C_M) / \ln 2$;
1.	10	0	$10 \cdot 0,5 / 12 = 0,417$	19	$1/19 = 0,0526$	sākotnējais $t_C^0 = 0$ s
2.	8	2	$8 \cdot 0,5 / 12 = 0,333$	24	$1/24 = 0,04167$	1,814062
3.	6	4	$6 \cdot 0,5 / 12 = 0,25$	33	$1/33 = 0,0303$	4,126011
4.	4	6	$4 \cdot 0,5 / 12 = 0,167$	55	$1/55 = 0,01818$	7,379777
5.	2	8	$2 \cdot 0,5 / 36 = 0,0833$	114	$1/114 = 0,00877$	12,98901

Write this result of reaction time t on table for test tubes (§ 1, § 2, § 3, § 4, § 5)!

Calculate $\vec{v} = 1/t \text{ sec}^{-1}$ reaction rate from data ((§ 1, § 2, § 3, § 4, § 5) in table!

Draw the graphic of velocity \vec{v} , sec^{-1} dependence on concentration $\text{CH}_2\text{S}_2\text{O}_3$

On graph calculate constant k in ratio \vec{v} , s^{-1} over concentration $\text{CH}_2\text{S}_2\text{O}_3 = 0,5 \text{ M}$:

Velocity $\vec{v} = k \cdot \text{CH}_2\text{S}_2\text{O}_3$; velocity constant: $k = \vec{v} / \text{CH}_2\text{S}_2\text{O}_3 = 0,062 / 0,5 = 0,124 \text{ s}^{-1}$

Half life time: $\tau^{1/2} = \frac{\ln 2}{k} = 0,693 / 0,124 = 5,59 \text{ s}$; $t_{1/10} = \frac{\ln 10}{k} = 2,3 / 0,124 = 18,6 \text{ s}$;


Secinājumi

a) Active mass law! Velocity \vec{v} is proportional to reactant concentration $\text{CH}_2\text{S}_2\text{O}_3$.

b) Velocity constant physical meaning $k = 0,124 \text{ s}^{-1}$ is velocity if concentration equal to $C = 1 \text{ M}$ one!

c) Half life time $\tau^{1/2} = 5,59$ seconds are during compounds amount decreases per half 1/2 !

d) Time $t_{1/10} = 18,6$ in seconds during thio sulfuric acid concentration $\text{CH}_2\text{S}_2\text{O}_3$ decreases 10 times?



Conclusions questions

- Active mass law! Velocity \vec{v} is proportional to reactant concentration $\text{CH}_2\text{S}_2\text{O}_3$.
- Velocity constant physical meaning $\vec{k} = 0,124 \text{ s}^{-1}$ is velocity if concentration equal to $C=1 \text{ M}$ one!
- Half life time $\tau_{1/2} = 5,59$ seconds are during compounds amount decreases per half $1/2$!
- Time $t_{1/10} = 18,6$ in seconds during thio sulfuric acid concentration $\text{CH}_2\text{S}_2\text{O}_3$ decreases 10 times?

Catalysis


Diverse catalysis in each of two test tubes pore approximately 1 mL H₂O₂ solution.

First test tube add in spoon of oxide Fe(+3)₂O₃. Secon test tube in add spoon of oxide Mn(+2)O.

1) Protolysis H₂O₂ deprotonate to anions 2H₂O₂+2H₂O→H₃O⁺+HOO⁻+OOH⁻+ H₃O⁺ pKa= 11,75.

2) Deprotonate negative anions HOO⁻><OOH⁻ collision energy E_a=79000 J/mol disproportion , oxidize colliding oxygen atoms to O₂aqua and second atoms reduce to 2 OH⁻ ions, which


3) neutralized to water 2OH⁻+2H₃O⁺=>2H₂O+2H₂O.



1) Reaction start by protolysis , than 2) collision disproportion OO atoms and 3) finally neutralized.



Negative ion collision with positive +2 HOO⁻>Mn(+2)O<OOH⁻ activation energy is E_a=790 J/mol



Catalase negative and positive ions collision HOO⁻>Fe³⁺ activation energy E_a=29 J/mol decreases:



Protolysis activate active transition state complex **oxygene+ water+ heat+ CATALYST**

Negative ions collision with positive iron Fe³⁺ ion create life resources 30 million times faster.

Describe observations and reaction equation. State, which catalyst effective disproportionate hydrogen peroxide H₂O₂. Determine, are given catalysts homogenous or heterogenic.

2. Autocatalyst (self-catalyst)

Each of two test tubes pore approximately 2 mL potassium permanganate KMnO₄ solution.

Both test tubes add approximately 1 mL sulfuric acid H₂SO₄ solution.








First test tube add spoon solid mangase(II) sulfate MnSO₄.

Both test tubes add approximately 1 mL oxalic acid H₂C₂O₄ solution.

Describe observations. Remark, in which test tube KMnO₄ collor disappears first.

Complete oxidation-reduction equation.

Izskaidrojiet mangāna(II) jonu Mn²⁺ lomū reakcijā!



Conclusions

1. Heterogeneous catalyst serves into reaction with surface area active centers.

are involved, therefore velocity of reaction is proportional to value surface area S.

2. Auto catalysis **KMnO₄ + H₂C₂O₄ + H₂SO₄** catalyses reaction product ion **Mn²⁺**.....